IVD

## 3diag - B2m - 800

REF TD-42521

LOT **521-16(1)**

2024-10

UDI

(01) 08434477102028  
(10) 521-16(1)  
(17) 241030  
(240) B28039

**EN** Beta-2 Microglobulin (serum+urine) - for *IMMAGE*<sup>®</sup> 800

**INSTRUCTIONS FOR USE**

**ES** Beta-2 Microglobulina (suero+orina) - para *IMMAGE*<sup>®</sup> 800

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**IT** Beta-2 Microglobulina (siero+urina) - per *IMMAGE*<sup>®</sup> 800

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**PT** Beta-2 Microglobulina (soro+urina) - para *IMMAGE*<sup>®</sup> 800

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<https://www.3diag.com/002>

for IFU, Scan or follow link,  
and select LOT

*IMMAGE*<sup>®</sup> is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

**TRIMERO Diagnostics, SL**

c. València 558, 4t 2a - 08026 Barcelona (Spain)  
☎ +34 93 244 86 79 - [www.3diag.com](http://www.3diag.com)

## INSTRUCTIONS FOR USE

Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)

# 3diag - B2m - 800

Beta-2 Microglobulin (serum+urine)

for **IMMAGE® 800**

**REF TD-42521**

(Product included in **REF TD-42520**)

### INTENDED USE

Quantitative determination of Beta-2 Microglobulin (B2m, B2μ), in human serum and urine, on **IMMAGE® 800** Immunochemical Systems of *Beckman Coulter* (**IMMAGE®** is a registered trademark of *Beckman Coulter, Inc.*, Fullerton, CA).

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- **UDR IMMAGE®** Reagent Cartridge: **REAG UDR B2m**  
**REF TD-42521-UDR** ▼ 100 test containing:
  - Antiserum Reagent: **CONT** 3.5 ml (Compartment A)  
Anti-human B2m antibodies bound to polystyrene particles.
  - Accelerator Reagent: **CONT** 5.5 ml (Compartment B)  
Solution of Ag/Ab reaction accelerators.

The reagents are ready for use and require no preparation. Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.

- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. When using the reagents, it is advisable to use anti-evaporation rubber stoppers (**REF** 447169, from *Beckman Coulter*). The indicated shelf life must be taken as a guideline given that, obviously, it depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- *Beckman Coulter's* **IMMAGE® 800** Immunochemical System, and accessories: dilution segments, cuvettes, etc..
- **BUF2 IMMAGE®**, from *Beckman Coulter* **REF** 447040
- **DIL1 IMMAGE®**, from *Beckman Coulter* **REF** 447640
- **3diag - B2m - Cal Set** **REF** TD-42522
- **3diag - B2m - Control (x3)** **REF** TD-42523
- **3diag - B2m - Control** **REF** TD-42533

### SAMPLES

- Fresh Serum.  
Samples with presence of fibrin should be centrifuged.  
Do not use hemolyzed, lipemic or contaminated samples.
- Fresh urine.  
Samples should be centrifuged.  
In bibliography<sup>(1)</sup> it is reported the following stability :
  - Serum: Refrigerated: 14 days
  - Urine: Refrigerated: 48 hours - Frozen: 60 daysThe pH must be adjusted between 6 and 8, since the B2M degrades rapidly if pH <6. It is recommended<sup>(1)</sup> to use 1M NaOH.

### PROCEDURE

Follow the instructions<sup>(2)</sup> for use of the analyzer to program and calibrate **UDR** (User Defined Reagent) assays, with the recommended parameters that are detailed below.

#### Assay Parameters - Serum/Urine

- Units: mg/l or IU/ml
- Protocol: Non-competitive NIPIA
- Sample or Dilution Volume: 30 μl
- Reaction Buffer Volume: 195 μl (use **BUF2**)
- Compartment A Volume: 30 μl
- Compartment B Volume: 50 μl
- Gain: 1
- Calibration Dilution: 1:30 (use **DIL1**)
- Sample Dilution: **Serum:** 1:30 (use **DIL1**)  
**Urine:** 1:5 (use **DIL1**)
- Reaction Time: 3.0 minutes

Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with Physiological Solution or **DIL1**, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

#### Calibration Parameters

- Calibrators: Use the **3diag - B2m - Cal Set**
- Levels: 6
- Replicates: 2
- Model: Logistic 4 Parameters

The assays must be recalibrated, at least when a new batch of reagents is used.

The use of the calibration update option with a single calibrator level is discouraged.

## PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assays are given in the Technical Reports, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

## QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of **3diag - B2m - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

## TRACEABILITY

Values are referred to the *1st International Standard for Beta2 Microglobulin* (code: *B2M*) of the *WHO (World Health Organisation)*.

## REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography reports the following reference values:

- Serum: between 1.1 and 2.4 mg/l<sup>(3)</sup>
- Urine: an excretion greater than 0.37 mg/24h<sup>(3)</sup> or a concentration of 0.3 mg/l can be considered significant<sup>(1)</sup>.

## CLINICAL SIGNIFICANCE

Elevated serum levels of B2m are either due to an increased production or to a decrease in the glomerular filtration rate which determines its accumulation. The increased rate of synthesis may be caused by many different diseases, in particular those of the lymphoid immune system, hematologic malignancies, viral infections, autoimmune diseases and others. B2m is hence not used as a specific diagnostic marker for any pathology in particular, but is rather widely used as a quantitative prognostic marker in many diseases. To take just one example, measurement of B2m levels is included in the International Staging System for Multiple Myeloma.

Measurement of B2m concentrations in urine is a very sensitive marker of tubular dysfunction due to any cause or pathology. Although as a tubular marker it suffers from the problems of instability in acidic environments (pH<6) and variable production levels.

## SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(4)</sup> by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

**REAG** Reagent

**UDR** UDR (User Defined Reagent)

**B2m** Beta-2 Microglobulin

**CONT** Contents

## BIBLIOGRAPHY

- (1) *Mayo Medical Laboratories* website ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), date of consultation: 11<sup>th</sup> July 2017.
- (2) "*IMMAGE® Immunochemical System user's manual*". *IMMAGE®* is a registered trademark of *Beckman Coulter, Inc.*, Fullerton, CA.
- (3) "*Nomenclator de Laboratorio Clinico AEFA/AEBM*" - ISBN: 84-486-0117-3.
- (4) EDMA Labelling Task Force: "*EDMA Symbols for IVD Reagents and Components - Revision, October 2009*".

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Os valores estão referidos ao *1st International Standard for Beta2 Microglobulin* (code: *B2M*) da *WHO (World Health Organisation)*.

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
**NOTE #1:** Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

**ALL PROTEINS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.