

3diag - B2m - 800

Beta-2 Microglobulin - for *IMMAGE*[®] 800

REF TD-42521

IVD CE

IFU Annex: *IMMAGE*[®] and *IMMAGE*[®] 800 Immunochemical Systems Application (Urine)

PROTOCOL INFORMATION

Reagent Lot Number:	(*1)	Units:	mg/L or IU/mL
Cartridge Lot Number:	(*1)	Protocol:	Non-competitive NIPIA
Cartridge Serial Number:	(*1)	AGXS Limit:	not defined
Reagent Expiry date:	(*1)		
Sample or Dilution Volume:	30 µL	Diluent:	DIL1 (*2)
Reaction Buffer Volume:	195 µL	Buffer:	BUF2 (*3)
Calibration Dilution:	1 : 30	Gain:	1
Compartment A Volume:	30 µL	Sample Dilution:	1 : 5 (*4)
Compartment B Volume:	50 µL	Reaction Time:	3.0 minutes

CALIBRATION INFORMATION

Levels:	6	Replicates:	2
Update Level:	4 (*5)	Replicates:	2 (*5)
Model:	Logistic 4 parameters		

NOTES

- (*1): See the labels of the Reagent Cartridge.
- (*2): Beckman Coulter P/N no. 447640. For the *IMMAGE*[®] it has to be defined as DIL10.
- (*3): Beckman Coulter P/N no. 447040. For the *IMMAGE*[®] it has to be defined as BUF10.
- (*4): It has to be programmed after performing the calibration process for the first time.
- (*5): The use of the calibration update, with only one calibrator level, is discouraged.