## 3diag - C1In - TIA

REF TD-42591

LOT **591-04(1)**

2024-11

UDI

(01) 08434477305030  
(10) 591-04(1)  
(17) 241130  
(240) TD-42591

**EN** C1 (Esterase) Inhibitor - for Turbidimetry  
**INSTRUCTIONS FOR USE**

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<https://www.3diag.com/005>

for IFU, Scan or follow link,  
and select LOT

**TRIMERO Diagnostics, SL**

c. València 558, 4t 2a - 08026 Barcelona (Spain)  
☎ +34 93 244 86 79 - [www.3diag.com](http://www.3diag.com)

## INSTRUCTIONS FOR USE

Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)

# 3diag - C1In - TIA

**C1 (Esterase) Inhibitor  
for Turbidimetry**

**REF TD-42591**

### INTENDED USE

Quantitative determination of Complement C1 (Esterase) Inhibitor (C1In, CEI), in human serum, by turbidimetric method in automatic Clinical Chemistry Analyzers.

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG Ab C1In**  
**REF** TD-42591-RA     ▽ 100 test - 6 ml  
Solution of anti-human C1In antibodies.
- Reaction Buffer: **BUF C1In**  
**REF** TD-42551-BF     ▽ 100 test - 18 ml  
TRIS Buffer, with PEG.

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- Automatic Clinical Chemistry Analyzer, capable of running photometric assays at 340 nm, and accessories: reagent containers, cuvettes, etc..

- 3diag - C1In - CAL SET**     **REF** TD-42582
- 3diag - C1In - CONTROL**     **REF** TD-42583

### SAMPLES

Fresh Serum.

Samples with presence of fibrin should be centrifuged.

Do not use hemolyzed, lipemic or contaminated samples.

In bibliography<sup>(1)</sup> it is reported the following stability in serum:

Refrigerated: 72 hours - Frozen: 14 days

### PROCEDURE

If necessary, carefully transfer the reagents to the containers used by the analyzer, preventing leakage and foaming or bubbles.

To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below. Please, contact the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79) for further information about applications to specific analyzers.

For some analyzers, in order to process the calibrators it may be necessary to deactivate the clot detection system.

### Assay Parameters

- Dispense and Mix:
  - Sample/Calibrator/Control: 20 µl (diluted 1:5)
  - BUF C1In** 180 µl
- Incubate a fixed time between 1 and 5 minutes
- Dispense and Mix:
  - REAG Ab C1In** 60 µl
- Read absorbance A1 (Blank) at 340 nm
- Incubate a fixed time of about 5 minutes
- Read absorbance A2 (End Point) at 340 nm
- Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators
- Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with Physiological Solution, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range

### Calibration Parameters

- Use the **3diag - C1In - CAL SET**.
- If the analyzer allows it, it is recommended to program two replicates of each calibration point.
- The calibrations are Non-linear. For the calculation it is recommended to use a 3<sup>rd</sup> Order Polynomial, a Logit or a Polygonal adjustment.
- The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

### PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

## QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of the **3diag - C1In - CONTROL**.

For some analyzers, in order to process the controls it may be necessary to deactivate the clot detection system.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

## TRACEABILITY

Given that certified reference materials are not available, the values are referred to internal standards of highly purified proteins. Traceability is ensured by measuring the C1In in the *International Ref. Preparation for human serum complement Factors (NIBSC code: W1032)* of the WHO (World Health Organization).


## REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography reports, depending on the method used and the population analyzed, reference values of between 19 to 37 mg/dl<sup>(1)</sup>, 15 to 35 mg/dl<sup>(2)</sup> and 17.4 to 24.0 mg/dl<sup>(3)</sup>.

Analyzing serum samples of 179 patients from Spanish hospital laboratories, without discarding those that presented alterations of other components of the complement, the following results have been obtained (in mg/dl):

media	DS	rango	percentil 95	percentil 90
20.7	3.91	9.18 - 38.5	14.2 - 28.5	15.0 - 26.8

In view of the results, a concentration lower than about 15 mg/dl can be taken as a significant value, indicative of a deficiency, since high C1In values do not have an established clinical significance<sup>(1)(3)</sup>.

## CLINICAL SIGNIFICANCE

C1 inhibitor deficiency is associated with hereditary or acquired angioedema. In 85% of cases their levels are decreased (5-30% of their normal value), while in the remaining 15%, the levels are normal but the protein is not functional. The additional measurement of C1q is key to distinguish between hereditary or acquired angioedema, since it presents normal levels in the hereditary form and low levels in the acquired form.

Low levels of C1 inhibitor also predispose to autoimmune diseases, especially Lupus Erythematosus (SLE), due to its effect of consumption of C3 and C4.

## SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(4)</sup> by the EDMA (*European Diagnostic Manufacturers Association*), whose meaning is detailed below.

**REAG** Reagent

**Ab** Antibody / Antiserum

**BUF** Buffer

**C1In** C1 (Esterase) Inhibitor

## BIBLIOGRAPHY

- (1) Mayo Medical Laboratories website ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)), date of consultation: 8<sup>th</sup> May 2018.
- (2) Putnam F.W. - The Plasma Protein - II Edition.
- (3) "Nomenclator de Laboratorio Clinico AEFA/AEBM" - ISBN: 84-486-0117-3.
- (4) EDMA Labelling Task Force: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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**NOTE #1:** Applications available for **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

IMMAGE® and AU® are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
 BN™ and ADVIA® are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
 Optilite® is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
 Alinity and Architect related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
 cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Roetkreuz, Switzerland

**CONTROLS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.