



c. València 558, 4t 2a - 08026 Barcelona +34 93 244 86 79 - info@3diag.com

> NIF: ES B66739772 Inscrita R.M. Barcelona Tomo 45329, Folio 165, Hoja 484362

3diag - C1In - TIA

REF TD-42591 - C1 (Esterase) Inhibitor - for Turbidimetry

ANNEX to IFU: cobas® 6000 (c501)

Proposal of Application - using 3diag - C1In - CAL SET (REF TD-42582)

cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Roctkreuz, Switzerland

<u>ANALYZE</u>						
Assay:	2 Point End	Time/Point:		oint:	[10] [37] [70] [0] [0]	
Wavelenght (2nd/Pri.):	[Cancel] [340]	Sample Vol. (Norm.):		rm.):	[3] [0] [0] Dil.Fact. = 1:1	
Sample Volume (Dec.):	[25] [3] [100]	Sample Volume (Inc.):		Inc.):	[0] [0] [0] - Not Defined	
Diluent:	[951 - Saline] [1]	R1 Volume:		ume:	[108] [0] [Inactive]	
R2 Volume:	[0] [0] [Inactive]	R3 Volume:		ume:	[36] [0] [Inactive]	
Linearity Limit:	[0] [0] [0]	Prozone Limit:		imit:	Not Defined	
Absorvance Limit:	[32000] [Increase]	Cell Detergent:		gent:	[Detergent 1]	
Stirring Level:	[2]					
Stirring Settin	Setting (UP - LOW-M1/M2/M3): [Stirring] - [Stirring] [Stirring] [Stirring]			g] [Sti	rring]	
<u>CALIBRATION</u>						
Calibration Type:	[RCM] or [Spline]	Point:		oint:	[6]	
Span - Weight:	[6] ^(*1) - [0]	Update Type:		Гуре:	[none] [0] [0]	
SD Limit:	[999]	Duplicate Limit (%/Abs.):		Abs.):	[99] [32000]	
Sensitivity Limit:	[-99999] [99999]	S1 Abs. Limit:			[-32000] [32000]	
<u>RANGE</u>						
Unit:	[mg/dl] or [mg/l]	Decimal Places:		aces:	[1] or [2] - User Defined	
Data Mode:	[Active]	Automatic Rerun:		erun:	User Defined	
Technical Limit - Low:	(*2)	Technical Limit - High:			(*3)	
OTHER (Calibrators)						
Std. Code, Rack and Conc. (#16):	(*4)					
Sample Volume (#16):	[3] [0] [0] [3] [0] [0]	[3] [0] [0]	[3] [0] [0]	[3] [0] [0] [3] [0] [0]	
Dilution Factor (#16):	1:1 1:1	1:1	1:1	1:1	1:1	
NOTES						
(*1) The use of the calibration update, with only one calibrator level, is discouraged.						

- (*2) User Defined (equal to lowest calibrator value (recommended)).
- (*3) User Defined (equal to highest calibrator value (recommended)).
- (*4) User Defined, for the concentrations see the values of the Calibrator Set on the table of assigned values of the IFU.
- (*5) The definition of these parameters is based on the *Roche*'s document "*cobas c pack MULTI*" (Doc. Code: 04929578001V1 Version: 2006-08, V 1 Español). For any doubt or further information, please refer to this document.