## 3diag - C1q - T1A

REF TD-42551

LOT **551-14**

 2024-10

UDI

(01) 08434477303036  
(10) 551-14  
(17) 241030  
(240) TD-42551

- EN** C1q Complement - for Turbidimetry  
**INSTRUCTIONS FOR USE**
- ES** Complemento C1q - para Turbidimetría  
**INSTRUCCIONES DE USO**
- FR** Complément C1q - pour Turbidimétrie  
**INSTRUCTIONS D'UTILISATION**
- IT** Complemento C1q - per Turbidimetria  
**ISTRUZIONI PER L'USO**
- NL** Complement C1q - voor Turbidimetrie  
**GEBRUIKSINSTRUCTIES**
- PT** Complemento C1q - para Turbidimetria  
**INSTRUÇÕES DE UTILIZAÇÃO**

<https://www.3diag.com/003>  
for IFU, Scan or follow link,  
and select LOT

**TRIMERO Diagnostics, SL**

c. València 558, 4t 2a - 08026 Barcelona (Spain)  
☎ +34 93 244 86 79 - [www.3diag.com](http://www.3diag.com)

## INSTRUCTIONS FOR USE

Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)

# 3diag - C1q - TIA

**C1q Complement  
for Turbidimetry**

**REF TD-42551**

(Product included in **REF TD-42555**)

### INTENDED USE

Quantitative determination of C1q Complement (C1q), in human serum, by turbidimetric method in automatic Clinical Chemistry Analyzers.

### PRINCIPLE OF THE METHOD

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

### CONTENTS - COMPOSITION - PREPARATION

- Antiserum Reagent: **REAG Ab C1q**  
**REF** TD-42551-RA ▽ 100 test <sup>(\*1)</sup> - 5 ml  
Anti-C1q antibodies, bound to polystyrene particles.
- Reaction Buffer: **BUF C1q**  
**REF** TD-42551-BF ▽ 100 test <sup>(\*1)</sup> - 20 ml  
TRIS Buffer, with PEG.

Note (\*1): with the recommended assay general parameters.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

### WARNINGS - PRECAUTIONS

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.
- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.

- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

### STORAGE - SHELF LIFE

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. This information should be taken as a guideline given that, obviously, the shelf life depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

### MATERIALS NEEDED, NOT SUPPLIED

- Automatic Clinical Chemistry Analyzer, capable of running photometric assays at 540...600 nm, and accessories: reagent containers, cuvettes, etc..
- 3diag - C1q - CAL SET** **REF** TD-42542
- 3diag - C1q - CAL** **REF** TD-42552
- 3diag - C1q - CONTROL** **REF** TD-42543

### SAMPLES

Fresh Serum.

Samples with presence of fibrin should be centrifuged.

Do not use hemolyzed, lipemic or contaminated samples.

In bibliography<sup>(1)</sup> it is reported the following stability in serum:

- Refrigerated: 10 days
- Frozen: 29 days.

### PROCEDURE

If necessary, carefully transfer the reagents to the containers used by the analyzer, preventing leakage and foaming or bubbles.

To program and calibrate assays, follow the instructions for use of the analyzer used, with the recommended general parameters that are detailed below.

Please, contact the Customer Support Service (☎ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 2448679) for further information about applications to specific analyzers.

### Assay Parameters

- Dispense and mix:
  - Sample/Calibrator/Control: 10 µl (diluted 1:30)
  - BUF C1q** 200 µl
- Incubate a fixed time between 1 and 5 minutes
- Dispense and mix:
  - REAG Ab C1q** 50 µl
- Read absorbance A1 (Blank) at 540...600 nm
- Incubate a fixed time of about 5 minutes
- Read absorbance A2 (End Point) at 540...600 nm
- Interpolate the absorbance increment (A2-A1) of the samples and controls in the curve obtained with the calibrators
- Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range. **It is recommended to use Saline as diluent.**

As an alternative, reagents can be mixed as first step, and the sample dispensed as starter.

### Calibration Parameters

- Use the **3diag - C1q - CAL SET** or, if using the **3diag - C1q - CAL** program in the analyzer or prepare the following dilutions: 1:1, 2:3, 1:2, 1:3, 1:6 and 1:12 (100, 66.7, 50, 33.3, 16.7 and 8.33 %).
- It is recommended to use **Saline** as diluent.
- If the analyzer allows it, it is recommended to program two replicates of each calibration point.

- The calibrations are Non-linear. For the calculation it is recommended to use a 3<sup>rd</sup> Order Polynomial, a Logit or a Polygonal adjustment.
- The assay must be recalibrated, at least when a new batch of reagents is used or when its parameterization is changed.

#### PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assay is given in the Technical Report, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 2448679).

#### QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of the **3diag - C1q - CONTROL**.

For some analyzers, in order to process the controls it may be necessary to deactivate the clot detection system.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

#### TRACEABILITY

Values in IU/ml are referred to the *International Ref. Preparation for human serum complement Factors* (NIBSC code: W1032) of the WHO (World Health Organization).

For the values in mg/dl, given that certified reference materials in these units are not available, the values are referred to internal standards of highly purified proteins. Traceability is ensured by measuring the C1q in the WHO standard.


#### REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography reports very variable reference values between publications, depending on the method used and the population analyzed: between 5.0 to 8.6 mg/dl<sup>(1)</sup>, 10 to 25 mg/dl<sup>(2)</sup>, 5.8 to 7.2 mg/dl<sup>(3)</sup>, 12 to 22 mg/dl<sup>(4)</sup>, and 17.5 to 48.5 mg/dl<sup>(5)</sup>.

Analyzing, by turbidimetric method, serum samples from 100 presumably healthy adult patients from the Barcelona area, the following results have been obtained (see table and histogram):

units	mean	SD	range	95 percentile
IU/ml	127	16.7	47.4 - 154	<b>100 - 147</b>
mg/dl	32.7	4.30	12.2 - 39.5	<b>25.7 - 37.9</b>

In view of the results, a concentration lower than about 100 IU/ml, equivalent to about 25 mg/dl, can be taken as a significant value, since high values of C1q have no established clinical significance<sup>(1)(3)</sup>.

#### CLINICAL SIGNIFICANCE

C1q component of the complement system is a glycoprotein with an approximate molecular weight of 400 kDa, composed of 18 peptide chains divided in 3 sub-units of 6. Together with the components C1r and C1s, it constitutes the C1 complex.

C1q recognizes and binds to Fc fragments of IgG and IgM immunoglobulins bound to antigen, thus activating the classical complement pathway cascade. The conformational changes produced by immunoglobulin binding enzymatically activate C1r and C1s which then continue the complement cascade.

In addition to activation of the complement classical pathway, another function of C1q is the elimination of immune complexes and apoptotic cells from the body.

C1q congenital deficiency is extremely rare (a few dozen cases in which the majority of patients suffered from Systemic Lupus (SLE)). Sometimes acquired deficiency is due to the presence of anti-C1q auto-antibodies.

Its deficiency has a significant effect on the host defense mechanisms and in the elimination of immune complexes. It is normally associated with a high incidence of autoimmune and infectious diseases.

Reduced levels of C1q are often found in acquired angioedema, due to the hyperactivation of the classical pathway of complement.

#### SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(6)</sup> by the EDMA (European Diagnostic Manufacturers Association), whose meaning is detailed below.

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
**NOTE #1:** Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

**ALL PROTEINS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.