

## 3diag - C1q - TIA

### C1q Complement - for Turbidimetry

REF TD-42551 IVD CE

## IFU ANNEX:

### AU400 - AU480 - AU640 Application

Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:		C1Q		<	>	Type:	Serum		Operation:	Yes	
Sample:	Volume	10.0	µL	Dilution	0	µL	Pre-Dilution Rate:	30			
Reagents:	R1 Volume	200	µL	Dilution	0	µL	Min OD	Max OD			
	R2 Volume	50	µL	Dilution	0	µL	L	-2.0	H	2.5	
Wavelength:	Pri.	570	nm	Sec.	None		Reagent OD limit:				
Method:	FIXED			First L	-2.0	First H	2.5				
Reaction slope:	+			Last L	-2.0	Last H	2.5				
Measuring Point 1:	First	12		Last	27		Dynamic Range:	L	**	H	**
Measuring Point 2:	First			Last			Correlation Factor:				
Linearity:			%				A	1.000	B	0.000	
No Lag Time:							On-board stability period:	#			

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:		C1Q		<	>	Type:	Serum			
Value/Flag:	#	Level L:	#	Level H:	#					
Normal Ranges:		Age L		Age H						
	Sex	Year	Month	Year	Month	L	H			
o 1.	#	#	#	#	#	#	#	#		
o 2.	#	#	#	#	#	#	#	#		
o 3.	#	#	#	#	#	#	#	#		
o 4.	#	#	#	#	#	#	#	#		
o 5.	#	#	#	#	#	#	#	#		
o 6.	#	#	#	#	#	#	#	#		
7. None Selected						#	#	#		
8. Out of Range		L	H		#	#	#	#		
Panic Value:	#		#	Unit:	mg/dL		Decimal places:	#		

Calibration Specific											
General		ISE									
Test Name:		C1Q		<	>	Type:	Serum				
Calibration Type:		6AB		Formula:	See Note (1)		Counts:	#			
Process:											
	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H						
Point 1:	1		**	-2.0	2.5						
Point 2:	2		**	-2.0	2.5						
Point 3:	3		**	-2.0	2.5						
Point 4:	4		**	-2.0	2.5						
Point 5:	5		**	-2.0	2.5						
Point 6:	6		**	-2.0	2.5						
Point 7:											
1-Point Cal. Point:				o With CONC-0							
MB Type Factor:		#		Calibration Stability Period:		#					

# User defined  
 \*\* Lot specific, see the analytical values table of the IFU of the Calibrators  
 (1) User defined POLYGONAL (recommended) or SPLINE