

3diag - C5 - TIA

ANNEX to IFU: *Architect c* - Application Proposal

REF TD-42571 - C5 Complement - for Turbidimetry

using **3diag - C5 - CAL SET** (**REF** TD-42562)

Architect and related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA

GENERAL

GENERAL

Assay	C5 (*1)	Assay Type	Photometric
Assay Availability	Enabled	Run Controls for R. by	Lot

REACTION DEFINITION

Reaction Mode	End Up	Primary Wavelength	340	Second. Wavelength	Not Defined
Main Read Time	32 - 33	Blank Read Time	18 - 19	Last Read Required	33
Sample Blank Type	Self	Absorbance Range	Not Defined	Color Correction	Not Defined

REAGENT / SAMPLE

Reagent	42571	R1 Reagent Volume	98	R2 Reagent Volume	42
Diluent Name	Saline	R1 Water Volume	45	R2 Water Volume	Not Defined
R1 Dispense Mode	Type 0	R2 Dispense Mode	Type 1		

<u>Dilution Name</u>	<u>Sample Volume</u>	<u>Dil. Sample Vol.</u>	<u>Diluent Volume</u>	<u>Water Volume</u>	<u>Dilution Factor</u>	<u>Default Dilution</u>
Std 1:1 (*1)	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)	Yes (Mark as Default)
D1 1:5 (*1)	25	3.5	100	Not Defined	5.00 (Informative)	No (Not Mark as Default)
3rd Dilution Not Defined						

VALIDITY CHECKS

Reaction Check Type	None	Read Time A Range	Not Defined	Calculation Limit	Not Defined
Minimum Absorv.	Not Defined	Read Time B Range	Not Defined	Rate Linearity %	Not Defined
Maximum Absorv. Variation	Not Defined				

CALIBRATION

CALIBRATION / CALIBRATORS

Calibration Method	Spline (Recommended)	Calibrator Set	C5 CAL (*1)	Replicates	2 (Recommended)
--------------------	-----------------------------	----------------	--------------------	------------	------------------------

NOTE: The Concentration of the calibrator level **C5 CAL 1**, used to made the blank, **must be set equal to zero**.

VOLUMES

<u>Cal Level</u>	<u>Sample Volume</u>	<u>Dil. Sample Vol.</u>	<u>Diluent Volume</u>	<u>Water Volume</u>	<u>Dilution Factor</u>
Blank: C5 CAL 1	1.5	3.5	300	Not Defined	201.00 (Informative)
Cal 1: C5 CAL 2	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)
Cal 2: C5 CAL 3	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)
Cal 3: C5 CAL 4	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)
Cal 4: C5 CAL 5	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)
Cal 5: C5 CAL 6	3.5	Not Defined	Not Defined	Not Defined	1.00 (Informative)

3diag - C5 - TIA

ANNEX to IFU: *Architect c* - Application Proposal

INTERVALS

<i>Full Interval (hours)</i>	Not Defined ^(*2)	<i>Adjust Interval (hours)</i>	Not Defined		
<i>Adjust Type</i>	None ^(*3)	<i>Adjust Level</i>	Not Defined	<i>Default Ordering Type</i>	Not Defined

VALIDITY CHECKS

<i>Blank Absorv. Range</i>	Not Defined	<i>Span</i>	Not Defined	<i>Span Absorv. Range</i>	Not Defined
<i>Expected Cal Factor</i>	Not Defined	<i>Exp. Cal. F. Toler. %</i>	Not Defined	<i>Maximum Curve Fit</i>	Not Defined

RESULTS

<i>Low Linearity</i>	User Defined ^(*4)	<i>High Linearity</i>	User Defined ^{(*4) (*5)}
<i>Gender and Age Spec. Ranges</i>	User Defined		
<i>Name, Range & Rev. Required</i>	User Defined		
<i>Result Units</i>	mg/dl	<i>Decimal Places</i>	2 (Recommended)
<i>Correlation Factor</i>	1.0000	<i>Intercept</i>	0.0000

NOTES

- (*1) Proposal, User defined field.
- (*2) Calibration curves have a limited validity, which depends on the particular conditions of use. We recommend to disable the automatic control of the calibration interval, and re-calibrate when:
 - a new lot of reagents is used,
 - established internal quality control procedures do not deliver the expected results, or
 - after performing maintenance operations on the analyzer.
- (*3) The use of the calibration adjustment, with only one or two calibrator levels, is discouraged.
- (*4) Linearity Limits can be left undefined, as for non-linear calibrations the analyzer automatically controls and flags samples with signals higher than the highest calibrator. If the user wants to use the Linearity Limits, we recommend to define it as:
 - Low Linearity Limit equal to 2.0 (fixed value), and
 - High Linearity Limit equal to Cal Set (REF: TD-42562) Level-6 value.
- (*5) High Linearity Limit should be adjusted to the new calibrator values whenever a new lot of Cal Set (REF: TD-42562) with different values is used.
- (*6) The number of test per kit can be optimized if 2 kits are transferred into a single container.
- (*7) It is recommended to retest samples higher than the upper limit of the assay range at 1:5 dilution.