

3diag - RBP - TIA

REF TD-42671 - Retinol Binding Protein - for Turbidimetry

ANNEX to IFU: *cobas*® 6000 (c501)

Proposal of Application - URINE - using 3diag - RBP - CAL SET (REF TD-42662)

cobas® and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

<u>ANALYZE</u>			
Assay:	2 Point End	Time/Point:	[10] [37] [70] [0] [0]
Wavelength (2nd/Pri.):	[Cancel] [600]	Sample Vol. (Norm.):	[20] [10] [80] Dil.Fact.=1:5
Sample Volume (Dec.):	[3] [10] [117] Dil.Fact.=1:40	Sample Volume (Inc.): (*1)	[10] [0] [0] Dil.Fact.=1:1
Diluent:	[951 - Saline] [1]	R1 Volume (Bottle "a" of pack MULTI):	[120] [0] [Inactive]
R2 Volume (Bottle "b" of pack MULTI):	[0] [0] [Inactive]	R3 Volume (Bottle "c" of pack MULTI):	[20] [0] [Inactive]
Linearity Limit/Prozone Limit:	Not Defined	Absorbance Limit:	[32000] [Increase]
Cell Detergent:	[Detergent 1]	Stirring Level:	[2]
Stirring Setting (UP - LOW-M1/M2/M3):	Leave pre-defined configuration (All Fields = Stirring)		
Cassette Type A (Reagent/No. Test/Vol. Use): (*2)	a: [R1] [50] [6.0] - b: [Cancel] [0] [0] - c: [R3] [50] [1.0]		
<u>CALIBRATION</u>			
Calibration Type:	[Spline]	Point:	[6]
Span - Weight:	[6] (*3) - [0]	Update Type:	[none] [0] [0]
SD Limit:	[999]	Duplicate Limit (%/Abs.):	[99] [32000]
Sensitivity Limit:	[-99999] [99999]	S1 Abs. Limit:	[-32000] [32000]
<u>RANGE</u>			
Unit:	[mg/dl] Proposal, User Defined	Decimal Places:	[4] Proposal, User Defined
Data Mode:	[Active]	Automatic Rerun:	User Defined
Technical Limit - Low:	(*4)	Technical Limit - High:	(*4)
<u>OTHER (Calibrators)</u>			
Std. Code, Rack and Conc. (#1..6):	(*5)		
Sample Volume (#1..6):	[20] [10] [80] [20] [10] [80] [20] [10] [80] [20] [10] [80] [20] [10] [80] [20] [10] [80]		
Dilution Factor (#1..6):	1:5	1:5	1:5
<u>NOTES</u>			
(*0)	The definition of these parameters is based on the Roche's document " <i>cobas c pack MULTI</i> " (Doc. Code: 04929578001V1 - Version: 2006-08, V 1 Español). For any doubt or further information, please refer to this document.		
(*1)	The definition of an Increased Sample Volume is optional, if not desired leave all fields = 0.		
(*2)	Calculated using the dead volumes and equations detailed in the Instructions of the " <i>cobas c pack MULTI</i> ". The number of test per kit could be optimized up to 165 test/pack if 2 kits are transferred into a single container.		
(*3)	The use of the calibration update, with only one calibrator level, is discouraged.		
(*4)	User Defined (Low Limit equal to Level-1 calibrator value and High Limit equal to Level-6 calibrator value - recommended).		
(*5)	User Defined, for the concentrations (lot specific) see the values of the 6-Level Calibrator Set (REF: TD-42662), on the table of assigned values of the IFU.		