## 3diag - sTfR - 800

REF TD-42681

LOT **681-13**

⏳ 2024-10

UDI

(01) 08434477110023  
(10) 681-13  
(17) 241030  
(240) B28074

**EN** Soluble Transferrin Receptor - for *IMMAGE*<sup>®</sup> 800  
**INSTRUCTIONS FOR USE**

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<https://www.3diag.com/010>  
for IFU, Scan or follow link,  
and select LOT

*IMMAGE*<sup>®</sup> is a registered trademark of Beckman Coulter, Inc., Fullerton, CA

**INSTRUCTIONS FOR USE**Reagents for professional use,  
for *In Vitro* use only in clinical laboratory (IVD)**3diag - sTfR - 800****Soluble Transferrin Receptor  
for IMMAGE® 800****REF TD-42681**

(Product included in REF TD-42680)

**INTENDED USE**

Quantitative determination of Soluble Transferrin Receptor (sTfR), in human serum, on IMMAGE® 800 Immunochemical Systems of Beckman Coulter (IMMAGE® is a registered trademark of Beckman Coulter, Inc., Fullerton, CA).

**PRINCIPLE OF THE METHOD**

The specific antibodies (Ab) of the reagent, bound to polystyrene particles, when combined with the antigens (Ag) of the patient sample, form insoluble compounds causing a change in the absorbance and dispersion of the light, proportional to the antigen concentration, which can be quantified by turbidimetric (TIA) or nephelometric (NIA) method, by comparison with calibrators of known concentration.

**CONTENTS - COMPOSITION - PREPARATION**

- UDR IMMAGE® Reagent Cartridge: **REAG UDR sTfR**  
REF TD-42681-UDR ▼ 100 test  
containing:
  - Antiserum Reagent: **CONT** 5.5 ml (Compartimento A)  
Solution of nnti-human sTfR antibodies, bound to polystyrene particles.
- Reaction Buffer: **BUF sTfR**  
REF TD-42681-BF ▼ 100 test - 35 ml  
TRIS Buffer, with PEG.

The reagents are ready for use and require no preparation.

Before each use it is convenient that the reagents are homogenized, shaking them gently avoiding the formation of foam or bubbles.

As a preservative, the reagents contain <0.1% (1 g/l) Sodium Azide (NaN<sub>3</sub>).

**WARNINGS - PRECAUTIONS**

- Sodium Azide is toxic. Even if sodium azide is not harmful at the concentration present in the reagents, take the necessary precautions to avoid accidental ingestion or contact with the eyes.
- Sodium Azide can react with lead or copper to give explosive compounds. For disposal it is recommended to rinse with plenty of running water to avoid accumulation in drains.

- Since the absence of infectious agents can not be proven with absolute certainty, components containing materials of human or animal origin must be handled with caution, as potentially infectious, following the recommended safety standards for biological risk.
- Do not mix components belonging to different lot kits.
- Clinical diagnosis should not be based on the results of a single test, but should always integrate all relevant clinical and laboratory data.

**STORAGE - SHELF LIFE**

- Store refrigerated at +2...+8°C. Do not freeze, as the functionality of the reagents may be altered.
- Properly stored and unopened, the reagents are stable until the expiration date indicated on the label.
- Once opened, the shelf life of the reagents is at least 4 weeks, provided that after each use they are stored immediately in the original containers, tightly capped and refrigerated at +2...+8°C. When using the reagents, it is advisable to use anti-evaporation rubber stoppers (REF 447169, from Beckman Coulter). The indicated shelf life must be taken as a guideline given that, obviously, it depends on the particular environmental and use conditions, which may differ from those of the stability studies carried out.

**MATERIALS NEEDED, NOT SUPPLIED**

- Beckman Coulter's IMMAGE® 800 Immunochemical System, and accessories: dilution segments, cuvettes, etc..
- DIL1 IMMAGE®, from Beckman Coulter REF 447640
- 3diag - sTfR - Cal Set REF TD-42682
- 3diag - sTfR - Control REF TD-42683

**SAMPLES**

Fresh Serum.

Samples with presence of fibrin should be centrifuged.

Do not use hemolyzed, lipemic or contaminated samples.

In bibliography<sup>(1)</sup> it is reported the following stability in serum:

- Refrigerated: 7 days
- Frozen: 90 days

**PROCEDURE**

Follow the instructions<sup>(2)</sup> for use of the analyzer to program and calibrate UDR (User Defined Reagent) assays, with the recommended parameters that are detailed below.

**Assay Parameters**

- Units: mg/l or nmol/l
- Protocol: Non-competitive NIPIA
- Sample or Dilution Volume: 10 µl
- Reaction Buffer Volume: 195 µl (use **BUF sTfR**)
- Compartment A Volume: 50 µl
- Compartment B Volume: 0 µl (not used)
- Gain: 1
- Calibration Dilution: 1:1 (neat)
- Sample Dilution: 1:5 (use DIL1)
- Reaction Time: 10.0 minutes

Samples with concentrations higher than the upper limit of the assay range should be analyzed again, diluted manually with Physiological Solution or DIL1, or by programming a larger sample dilution in the analyzer, to recover a value close to the midpoint of the measurement range.

**Calibration Parameters**

- Calibrators: Use the **3diag - sTfR - Cal Set**
- Levels: 6
- Replicates: 2
- Model: Logistic 4 Parameters

The assays must be recalibrated, at least when a new batch of reagents is used.

The use of the calibration update option with a single calibrator level is discouraged.

## PERFORMANCES OF THE METHOD

Detailed information on the characteristics and performances of the assays are given in the Technical Reports, available on the website ([www.3diag.com](http://www.3diag.com)) or upon request to the Customer Support Service (✉ [support@3diag.com](mailto:support@3diag.com) - ☎ +34 93 244 86 79).

## QUALITY CONTROL

To monitor performances, it is recommended that internal controls be inserted into each analytical series. It is recommended to use the controls of **3diag - sTfR - CONTROL**.

Each laboratory should establish its own quality scheme and corrective actions if controls do not meet the assigned tolerances.

The reagents have been subjected to quality control checks and should react as described in these instructions. Therefore, as a general recommendation, in case the controls do not give the expected reaction, as a precaution all reagents should be considered unreliable until their operation has been checked.

## TRACEABILITY

Values are referred to the *Reference Reagent Recombinant Soluble Transferrin Receptor (rsTfR)* (code: 07/202) of the *WHO (World Health Organization)*, and apply to free rsTfR monomer<sup>(3)</sup>.

## REFERENCE INTERVALS

It is always advisable for each laboratory to establish its own reference values.

The bibliography<sup>(4)</sup> reports reference values between 10.6 and 24.6 nmol/l, equivalent to 0.76 and 1.76 mg/l.

## CLINICAL SIGNIFICANCE

sTfR is a good marker in the diagnosis of iron deficiency. sTfR enables evaluation of erythropoiesis without the need for cytological study of the bone marrow, as long as iron deficiency is excluded. Thus, it is helpful in managing response to erythropoietin treatment (EPO), and it has also been proposed for use in anti-doping control.

The concentration of sTfR increases significantly and at an early stage in patients suffering from iron deficiency. Because it does not act as an acute phase reactant, it is especially helpful in the differential diagnosis of iron-deficiency anaemia and secondary anaemia from chronic disorders (not caused by cancer) in the presence of acute or inflammatory conditions (commonly seen in elderly patients for example) which affect Ferritin (Ft) measurements. Reports suggest that the sTfR/Log(Ft) ratio results in higher sensitivity and specificity of differentiation with respect to the individual magnitudes.

In its compilation of recommendations<sup>(5)</sup>, the *World Health Organization (WHO)* concludes that sTfR levels, combined with measurements of serum ferritin concentrations, provide the most effective assessment of a population's iron status.

## SYMBOLS

In addition to the harmonized symbols provided on the European Standard EN 980:2008, in the labels and instructions of use has been used the complementary symbology proposed<sup>(6)</sup> by the *EDMA (European Diagnostic Manufacturers Association)*, whose meaning is detailed below.

**REAG** Reagent

**UDR** UDR (User Defined Reagent)

**STfR** Soluble Transferrin Receptor

**CONT** Contents

**BUF** Buffer

## BIBLIOGRAPHY

- (1) *Mayo Medical Laboratories* ([www.mayomedicallaboratories.com](http://www.mayomedicallaboratories.com)) website, date of consultation: 30<sup>th</sup> May 2017.
- (2) *IMMAGE*® Immunochemical System user's manual. *IMMAGE*® is a registered trademark of *Beckman Coulter, Inc.*, Fullerton, CA.
- (3) For further information visit *NIBSC (National Institute for Biological Standards and Control)* website: [www.nibsc.org](http://www.nibsc.org).
- (4) *Quest Diagnostics*™ website ([www.questdiagnostics.com](http://www.questdiagnostics.com)), date of consultation: 17<sup>th</sup> June 2017.
- (5) "Serum transferrin receptor levels for the assessment of iron status and iron deficiency in populations. Vitamin and Mineral Nutrition Information System." Geneva: *World Health Organization*; 2014 (WHO/NMH/NHD/MNM/14.6; [http://apps.who.int/iris/bitstream/10665/133707/1/WHO\\_NMH\\_NHD\\_EPG\\_14.6\\_eng.pdf?ua=1](http://apps.who.int/iris/bitstream/10665/133707/1/WHO_NMH_NHD_EPG_14.6_eng.pdf?ua=1)), date of consultation: 30<sup>th</sup> September 2019.
- (6) *EDMA Labelling Task Force*: "EDMA Symbols for IVD Reagents and Components - Revision, October 2009".

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
**NOTE #1:** Applications available for **Optilite®**, **Alinity c**, **Architect c**, **ADVIA®** series, **AU®** series, **cobas®** series, and other analyzers by request

**IMAGE®** and **AU®** are registered trademarks of Beckman Coulter, Inc, Fullerton, CA.  
**BN™** and **ADVIA®** are registered trademarks of Siemens Healthcare Diagnostics Products GmbH, Marburg, Deutschland.  
**Optilite®** is a registered trademark of The Binding Site Group Ltd, Birmingham, U.K.  
**Alinity** and **Architect** related brand marks are registered trademarks of Abbott Laboratories, Abbott Park, Illinois, USA  
**cobas®** and related brand marks are trademarks of Roche Diagnostics Ltd, Rottkreuz, Switzerland

**ALL PROTEINS INCLUDED IN**  **Accuracy 365** Online Quality Control Further information and registration at [www.accuracy365.com](http://www.accuracy365.com)

**Accuracy 365** aims to monitor and compare, between laboratories, the results obtained in the Internal Quality Control processes of end-users of specific protein control materials manufactured by **TRIMERO Diagnostics**.

### Main features:

- 100 % free and lifetime.
- Cloud application, Available 24 hours, 365 days per year.
- Intuitive interface and very simple to use.
- Calculation and representation of statistical data in real time.
- Anonymous participation, data are only accessible by its owner.